

	SOP-BCR-5.17	Double Immunofluorescence Staining for RAD-51 and H2Ax Foci in Cells	Author: S. Clouthier  Approved: M. Wicha 	Rev: 0	Issued: 6/19/2013 Revised:
--	--------------	---	--	--------	-------------------------------

1.0 Purpose




The purpose of SOP 5.17 is to state the procedure double Immunofluorescence Staining for RAD-51 and H2Ax Foci in Cells.

2.0 Scope

SOP is intended to cover all resources, personnel and equipment in the BCR laboratory.

3.0 Procedure

- 3.1 Plate 200k-300k cells into each well of two-well glass chamber slides (Collagen-coated for SUM159)
- 3.2 Irradiate cells with γ -irradiation at 4 Gray.
- 3.3 Change media immediately after irradiation and return cells to incubator.
- 3.4 After 8-24 hours (HCC1954: 15 hours, SUM149: 8 hours, SUM159: 12 hours), wash cells with room-temperature 1x PBS. Repeat this step three times.
- 3.5 Fix cells in room-temperature (RT) 2% ParaFormaldehyde (diluted in PBS). Incubate at room temperature for 15 minutes.
- 3.6 Wash cells with room temperature 1x PBS. Repeat this step three times.
 - If necessary, cells can be stored in PBS at 4°C overnight.
- 3.7 Permealize in 0.1% Triton X-100 (in PBS) for five minutes at room temperature. (Can use small glass slide holder).
- 3.8 Wash cells with room temperature 1x PBS. Repeat this step three times.
- 3.9 Use the blue blocking solution from AEC Staining Kit (Invitrogen 85-8943) or 10% Goat Serum from Invitrogen. This step takes 10 minutes.
- 3.10 **H2AX Primary Antibody** (-20°C Ab box in freezer 1; Cat#16-193, Millipore JBW301)
 - 1:1,000 dilution in 1x PBS
 - Incubate for 1 hour at room temperature
- 3.11 Wash cells with room temperature 1x PBS. Repeat this step three times.
- 3.12 **Turn off the lights.** Proceed with the following steps in the **DARK**.
- 3.13 **Alexa-fluor 546 Goat Anti-Mouse (Red)** (Invitrogen A11030 located in Freezer #2)
 - 1:2,000 dilution in 1x PBS.
 - Incubate for 30 minutes.
- 3.14 Wash cells with room temperature 1x PBS. Repeat this step three times.
- 3.15 Use the **blue blocking solution** from AEC Staining Kit (Invitrogen 85-8943) or 10% Goat Serum from Invitrogen. This step takes 20 minutes.
- 3.16 **RAD-51 Primary Antibody** (4°C; SantaCruz sc-8349; located in Unconjugated Ab Box in Fridge #1)
 - 1:200 dilution in 1x PBS.
 - Incubate for 1 hour at room temperature.
- 3.17 Wash cells with room temperature 1x PBS. Repeat this step three times.
- 3.18 **Alexa-fluor 488 Goat Anti-Rabbit (Green)** (Invitrogen A11008; located in 2^o Ab Box in Fridge #1)
 - 1:1,000 dilution in 1x PBS (1:200 too much background).
 - Incubate for 30 minutes.
- 3.19 Wash cells with room temperature 1x PBS. Repeat this step three times.

	SOP-BCR-5.17	Double Immunofluorescence Staining for RAD-51 and H2Ax Foci in Cells	Author: S. Clouthier  Approved: M. Wicha 	Rev: 0	Issued: 6/19/2013 Revised:
--	--------------	---	--	--------	-------------------------------

3.20 **Mount slides with Invitrogen ProLong Gold Antifade Reagent with DAPI** (cat # P-36931; located in Freezer #4)

3.21 Store slides overnight in the **DARK** at room temperature.

3.22 Analyze slides and capture photos after overnight storage

4.0 Applicable References

5.0 Change Description

Revision	Date	Reference	Description of Change