

	SOP-BCR-5.7	<b>Staining procedure for Tissue Sections (Goat Primary Antibody) by using RTU Vectastain Universal Quick Kit</b>	Author: S. Clouthier  Approved: M. Wicha 	Rev: 1.0	Issued: 09/24/98 Revised: 7/17/12
--	-------------	---	--	----------	--------------------------------------

## 1.0 Purpose

The purpose of SOP 5.7 is to provide instructions on how to stain for tissue sections (Goat primary antibody) by using RTU Vectastain Universal Quick Kit.

## 2.0 Scope

SOP 5.7 is intended to cover all resources, personnel and equipment in the BCR laboratory.

## 3.0 Materials

No.	Name	Description	Storage Location
1.0	1x PBS	Phosphate buffering saline	Cold Storage (026-380C)
2.0	Peroxidase quenching solution	(3% hydrogen peroxidase in methanol)	
3.0	RTU Biotinylated pan specific universal secondary antibody	Binds to primary antibody	Fridge #1 (026-328S-A)
4.0	Streptavidin/oxidase complex reagent	Binds to any biotinylated antibody	Fridge #1 (026-328S-A)
5.0	Peroxidase substrate solution	Staining Agent	Freezer #2 (026-328S-A)

## 4.0 Procedure

4.1 Deparaffine the slides cut from paraffin-embedded Tissue section.

- 3x xylenes for > 15 minutes
- 2x 100% EtOH for 5 minutes each
- Wash in dH<sub>2</sub>O for > 5 minutes in running water.

4.2 Antigen unmasking (Retrieve): Soak the slides in ready to use Citrate Buffer (pH 6.0) and heated at 98° C for 40 minutes and let it cool down for about 10 minutes (to room temperature). Wash in running dH<sub>2</sub>O for > 15 minutes.

4.3 Submerge slides in peroxidase quenching solution (3% hydrogen peroxide in methanol) for 10 minutes.

4.4 Wash slides in 1x PBS for 5 minutes

4.5 Incubate the slides in blocking serum for 10 minutes.

4.6 Blot excess serum from sections.

4.7 Incubate sections in primary antibody diluted in PBS containing 1.5% blocking serum for 1 hour (For **Scribble: 1:100**)




4.8 Wash slides for 5 minutes in PBS.

4.9 Incubate sections in RTU biotinylated pan specific universal secondary antibody for 10 minutes.

4.10 Wash sections for 5 minutes with PBS.

4.11 Incubate sections in ready-to-use streptavidin/oxidase complex reagent for 5 minutes.

4.12 Wash sections for 5 minutes with PBS.

	SOP-BCR-5.7	<b>Staining procedure for Tissue Sections (Goat Primary Antibody) by using RTU Vectastain Universal Quick Kit</b>	Author: S. Clouthier  Approved: M. Wicha 	Rev: 1.0	Issued: 09/24/98 Revised: 7/17/12
--	-------------	---	--	-------------	--------------------------------------

- 4.13 Protect from the light!
- 4.14 Incubate sections in peroxidase substrate solution (DAB Substrate Kit for peroxidase) until desired stain intensity develops. Immediately before use, prepare the substrate solution as follows.
- To 5.0 mL of distilled water, add 2 drops of Buffer Stock Solution and mix well.
  - Add 4 drops of DAB Stock Solution and mix well.
  - Add 2 drops of the Hydrogen Peroxide Solution and mix well.
  - If a gray-black stain is desired, add 2 drops of the Nickel Solution and mix well.
  - Incubate tissue sections with the substrate at room temperature until suitable staining develops. Development times should be determined by the investigator but generally 2-10 minutes provides good staining intensity.
- 4.15 Wash sections in dH<sub>2</sub>O for 5 minutes.
- 4.16 Counter stain, clear and mount. (Optional).

## 5.0 Applicable References

## 6.0 Change Description

Revision	Date	Reference	Description of Change
1.0	7/17/12	CL	Added room locations