

How to Use Your Flovent® (fluticasone) Diskus Inhaler

50 mcg 100mcg 250 mcg



What is my prescription?	Special Instructions
	Take every day to control your symptoms.
	Rinse mouth thoroughly with water and spit after each use.

What can I expect while I take this medication?

Flovent® Diskus is an inhaled corticosteroid that decreases airway swelling. **Flovent® Diskus** will not provide immediate relief, but you can expect fewer episodes of coughing and wheezing.

How do I use my Flovent® Diskus?

- 1. Hold the diskus in one hand & put the thumb of your other hand on the notch & push towards the back until it snaps & the mouthpiece and lever are revealed.
- 2. Slide the lever away from you until it clicks. This loads the medication.
- 3. Turn your head away from the diskus and breathe out as much air as you comfortably can.
- 4. Put the mouthpiece in your mouth. Make sure your teeth and tongue do not block the mouthpiece.
- 5. Breathe in as steadily and deeply as you can. Remove the diskus from

your mouth while holding your breath for 5-10 seconds. Breathe out.

- 6. Close the diskus by sliding the thumb grip towards you
- 7. Rinse mouth thoroughly with water after taking medication to decrease the risk of a mouth infection. Spit out the water. **Do not swallow it.**
- 8. You may not taste, feel or smell the dry powder when you are using your diskus. As long as you are following the directions, you will get your full dose of medication.

Do not use Flovent® Diskus to treat sudden symptoms of your health condition. Always carry a rescue inhaler with you.

To learn more about your health condition visit: http://www.med.umich.edu/careguides.

To learn more about the **Flovent® Diskus** inhaler visit:

http://www.flovent.com/use/flovent-diskus.html

Disclaimer: This document contains information and/or instructional materials developed by the University of Michigan Health System (UMHS) for the typical patient with your condition. It may include links to online content that was not created by UMHS and for which UMHS does not assume responsibility. It does not replace medical advice from your health care provider because your experience may differ from that of the typical patient. Talk to your health care provider if you have any questions about this document, your condition or your treatment plan.

Patient Education by <u>University of Michigan Health System</u> is licensed under a <u>Creative Commons</u>
<u>Attribution-NonCommercial-ShareAlike 3.0 Unported License</u>. Last
Revised 4/15/2015

Flovent®Image: Copyright GlaxoSmithKline. Used with permission.

COPD QI Steering Committee How to use your Flovent® Diskus