



University of Michigan
C.S. Mott Children's Hospital

Our Recommendations for Blood Sugar Printouts

Useful for Review by Diabetes Team

Accu-Check 360®

- Trend Graph
- Blood Glucose Levels

Animas ezManager Max®

- Logbook with BG info, carbohydrate info, bolus info (Show All Values)
- Pump Settings

Bayer Glucofacts®

- In Setting Preferences - Change the number of time periods to 5 (overnight, breakfast, lunch, dinner, bedtime) and adjust the times accordingly.
- Logbook
- Standard Day

DexCom®

- Glucose Trend
- Hourly Statistics
- Glucose Distribution

CoPilot®

- Logbook
- Weekly Pump
- Statistics
- Glucose Modal Day

IQHealth®

- Diary Summary - Standard Day - All Meals - View Details

Medtronic CareLink®

- Quick View Summary
- Logbook Diary
- Modal Day Hourly
- Device Settings
- Sensor Daily Overlay (for CGM)
- Sensor Overlay by Meal (for CGM)

OneTouch® Diabetes Software

- Logbook
- Standard Day

Disclaimer: This document is for informational purposes only and is not intended to take the place of the care and attention of your personal physician or other professional medical services. Talk with your doctor if you have questions about individual health concerns or specific treatment options.

©2011 The Regents of the University of Michigan
Last Revised 01/2012