

# Medical Shipping Symbols



Batch code.



Do not reuse



Outer Diameter



Inner Diameter



Manufacturer



Date of Manufacturer



Sterilized Using Irradiation



Consult Instructions For Use

**F**

French Size



Use By



Corrosive (Acid)



Attention, See Instructions for Use



Biological Risks



Content (Numeral represents quantity of units inside)



Caution radioactive Material



For use within temperature limits.



For use above a minimum limit.



For use below a minimum limit.



Control Symbol.



Control Positive.



Control Negative.



Idem: FDA version from IDO 15223.



Idem: Showing sterilization method: Ethylene Oxide.



Idem: Aseptic Processing.




Idem: Radiation.



Idem: Heat.



For In Vitro Devices Performance Evaluation Only.



Authorized EC Representative in the European Community.



Store in this position only. Keep upright.



In Vitro Diagnostic Medical Device.



Sufficient for N tests.




Caution, consult accompanying documents.

SN

Serial Number



Contents (Numeral represents quantity of units inside.)



Stent Length

**REF**

Catalog or reference number.



Ce mark