

Nikon Diaphot 200 Fluorescent Filter Cubes

Alternate 25mm Tray

Slider position	Use	Excitation	Dichroic Mirror	Emission
RIGHT (I.D. UV-2A)	UV	330 - 380 nm	400 nm	420+ nm LP
CENTER (I.D. HQ EN GFP)	eGFP	450 - 490 nm	495 nm	500 - 550 nm
LEFT (I.D. Cy5T)	Cy-5 (far red)	620 - 660 nm	660 nm	665 - 690 nm