

Nikon Diaphot 200 Fluorescent Filter Cubes

Tray 1

Slider position	Use	Excitation	Dichroic Mirror	Emission
RIGHT (I.D. UV-2E/C)	DAPI/Hoechst	340 - 380 nm	400 nm	435 - 485 nm
CENTER (I.D. B-2E/C)	FITC (also eGFP)	465 - 495 nm	505 nm	515 - 555 nm
LEFT (I.D. Y-2E/C TR)	Texas Red (also Rhodamine)	540 - 580 nm	595 nm	600 - 660 nm