



# HeartMate 3™ Patient Daily Log

Date							
Pump Speed (rpm)							
Flow Rate (lpm)							
Pulsatility Index (PI)							
Power (Watts)							
Check when daily controller self-test is done							
Describe the driveline site: (drainage, warmth, redness, swelling, odor)							
Blood Pressure							
Temperature							
Weight							

## Glossary of Terms:

---

Heartmate™ Term	Definition
<b>Pump Speed</b>	How fast the LVAD pump is spinning to circulate blood through your body; it is measured in revolutions per minute (RPM).
<b>Flow Rate</b>	The amount of blood the LVAD is circulating through your body; it is measured in liters per minute (LPM).
<b>Pulsatility Index (PI)</b>	The PI value is related to the amount of assistance provided by the blood pump
<b>Power</b>	The amount of energy it takes to power your pump at the speed at which it is set; it is measured in watts.

Disclaimer: This document contains information and/or instructional materials developed by Michigan Medicine for the typical patient with your condition. It may include links to online content that was not created by Michigan Medicine and for which Michigan Medicine does not assume responsibility. It does not replace medical advice from your health care provider because your experience may differ from that of the typical patient. Talk to your health care provider if you have any questions about this document, your condition or your treatment plan.

Reviewers: Bethany Lee-Lehner, RN, MSN  
CVC #1106

Patient Education by [Michigan Medicine](#) is licensed under a [Creative Commons Attribution-NonCommercial-ShareAlike 3.0 Unported License](#). Last Revised 3/2020