

	SOP-BCR-4.14	Preparation of Fungizone	Author: S. Clouthier  Approved: M. Wicha 	Rev: 1.0	Issued: 09/24/98 Revised: 7/17/12
--	--------------	---------------------------------	--	----------	--------------------------------------

1.0 Purpose

The purpose of SOP 4.14 is to describe how to prepare Fungizone

2.0 Scope

SOP 4.14 is intended to cover all resources, personnel and equipment in the BCR laboratory.

3.0 Materials & Reagents

No.	Name	Description	Storage Location
1.0	Fungizone	Amphotericin B	Freezer #2 (026-328S-A)

4.0 Procedure

- 4.1 Source: Cellgro 30-0003-CF
- 4.2 Thaw in 37°C water bath
- 4.3 Add 5 mL to 500 mL medium
- 4.4 Store at **4°C** after opening, store at **-20°C** before use

5.0 Applicable References

6.0 Change Description

Revision	Date	Reference	Description of Change
1.0	7/17/12	CL	Added room location